

**KAREI™****PTFE Series**

## DESCRIPTIONS

**PTFE** membrane cartridge filters are made of GORETEX® expanded Polytetrafluoro-Ethylene (ePTFE) media manufactured by DUPONT and supported with 100% pure Polypropylene material parts.

Inert Teflon and PP construction provides excellent chemical compatibility.

Expanded hydrophobic with high void volume properties of **PTFE** provides superior flow rate, low pressure drops and great retention capability. Make it the ideal choice for broad range of liquid and bulk air filtrations.

All parts are thermally welded without surfactants, additives and binders eliminate extractable.

3-Dimensional bonded continuous fibrils provide fixed pore structure. Eliminates shedding and particles unloading.

Excellent retention capability of micro-fine particles and bacteria (In specific pore sizes).

100% integrity tested to ensure product consistency.

**PTFE** is pre-rinsed with 18MΩ-cm D.I. water (Optional).

Available in a wide range of absolute pore size of 0.05, 0.1, 0.2, 0.45, 0.7 and 1 micron

Manufacture in Class 10K clean room environment minimizes contamination.

Comply with FDA Code Of Federal Regulation Title 21 for food and beverage use.

Meet USP Class VI-121°C Plastic reactivity test for Biosafety.

Removal rating of >99.9999% at 0.03 um particles.

A guaranteed quality product (ISO 9001 certified).

Cartridges will be rinsed-up to 18 MΩ-cm D.I. water with a

## SPECIFICATIONS

### ABSOLUTE MICRON RATING

0.05, 0.1, 0.2, 0.45, 0.7 and 1.0 micron

### FILTRATION AREA/ 10" FILTER CARTRIDGE

0.05, 0.1 & 0.2 um : 1 m<sup>2</sup> (10.8 ft<sup>2</sup>)

0.45 um : 0.85 m<sup>2</sup> (9.2 ft<sup>2</sup>)

0.7 & 1.0 um : 0.7 m<sup>2</sup> (7.6 ft<sup>2</sup>)

Note: For double layer of 0.05, 0.1, 0.2 & 0.45 micron, filtration area : 0.7 m<sup>2</sup> (7.6 ft<sup>2</sup>)

### NOMINAL LENGTH

125, 250, 500, 750, 1000 mm or

127, 254, 508, 762, 1016 mm

### NOMINAL INNER/OUTER DIAMETER (ID/OD)

Standard : 30/ 68 mm or BB : 28/ 114 mm

Note: 28mm inner diameter is available upon request.

### MEDIA MATERIAL

- i) HB : Hydrophobic Expanded PTFE membrane (For gases applications)
- ii) HL : Hydrophilic Expanded PTFE membrane (For liquid applications)

### SUPPORTING MATERIAL

Pure Polypropylene & Polyethylene Microfiber or PTFE

### INNER CORE, CAGE AND END ADAPTOR MATERIAL

- i) Standard : High Strength Pure Polypropylene
- ii) RPG : Reinforced Polypropylene With Glass
- iii) HPE : High Density Polyethylene
- iv) PTFE : PTFE

### SEALING TECHNIQUE

Thermal Bonding

### END STYLE

- 1) DOE : Double Opened End
- 2) SOE : Single Opened End
- i) S2C : SOE, 222 O-Ring With Closed End
- ii) S2F : SOE, 222 O-Ring With Finned End
- iii) S6C : SOE, 226 O-Ring With Closed End
- iv) S6F : SOE, 226 O-Ring With Finned End

Note: Extended adaptor for SOE filter cartridge is available upon request.

### GASKET AND O-RING MATERIAL

- 1) Standard : EPDM
- 2) V : Viton
- 3) S : Silicone
- 4) T : Teflon
- 5) FEP : Teflon Encapsulated Viton

## OPERATING CONDITIONS

### MAX. FORWARD DIFFERENTIAL PRESSURE

2.1 Bar (30 PSI) At 95°C, 4.8 Bar (70 PSI) At 25°C

### MAX. REVERSE DIFFERENTIAL PRESSURE



3.4 Bar (50 PSI) At 25°C

### MAX. OPERATING TEMPERATURE

90°C at 2.1 Bar (30 PSI)

### CHANGE OUT DIFFERENTIAL PRESSURE

2.4 Bar (35 PSID)



Air Temperature: 24 °C, Inlet Pressure: 15 PSID (1 Bar)

## STERILIZATION AND SANITIZATION METHODS

<b>Chemicals</b>	Peracetic acid, chlorinated alkaline products, bleach, sulfur dioxide (SO <sub>2</sub> ) and hydrogen peroxide at typical sanitization temperatures and concentrations.
<b>Autoclave</b>	125 °C (257°F) for 30-45 minutes at maximum differential pressure of 7 PSI (0.5 Bar).
<b>In-line Stream</b>	140 °C (284°F) for 45-60 minutes at 2 PSID (0.14 Bar) ΔP.
<b>Hot Water</b>	88 °C (190°F) at 5 PSI (0.3 Bar) up to 50 minutes.

**KAREI-PTFE** can be subjected to multiple sterilization cycles while maintain its integrity.

## TYPICAL APPLICATIONS

Chemicals Industries	Process Gases	Food & Beverage Industries	Pharmaceutical Industries
Rinses And Solvents, Concentrated Acids & Bases, Ozonated D.I. Water, Etchants, Photo Resists, Wet-Etch System, Developers, Strippers, Polymer Filtration.	Bulk And Point-Of-Use Gases, Compressed Air	Sterile Venting Of Holding Tanks, Sterile CO <sub>2</sub> Filtration, Microbial Control Of Inlet Air For Bio-Processing Of Foods, Water Process With UV, etc.	Tank Vents, Filtration Of Compressed Gases, Filtration Of Solvents, Extraction Of Chemicals, Re-Crystallization Of Chemicals.

## INTEGRITY TEST – MINIMUM BUBBLE POINT (IPA)

Micron	0.05	0.1	0.2	0.45	0.7	1.2
<b>BAR</b>	>2.4	>1.7	>1.4	>0.6	>0.4	>0.3
<b>PSIG</b>	>35	>25	>20	>9	>6	>5.5
<b>Log Retention Value Of Bacteria</b>	-	-	≥7	≥7	≥7	-
<b>Bacteria</b>	-	-	Brevundimonas Diminuta	Serratia Marcescens	Sacch. Cerevisiae	-

# CHEMICAL COMPATIBILITY GUIDE

<b>Acids</b>	Acetic Acid, Glacial	R	<b>Glycol</b>	Ethylene Glycol	R	
	Acetic Acid, 90%	R		Glycerine	R	
	Acetic Acid, 30%	R		Propylene Glycol	R	
	Boric Acid	R		<b>Hydrocarbons</b>	Benzene	LR
	Hydrochloric Acid, Conc.	R			Toluene	LR
	Hydrochloric Acid, 6N	R			Xylene	LR
	Hydrofluoric Acid, 6N	R			Carbon Tetrachloride	LR
	Nitric Acid, Conc.	R			Chloroform	LR
	Phosphoric Acid, Conc.	R			Freon TF	R
	Sulfuric Acid, Conc.	R			Methylene Chlorine	LR
Sulfuric Acid, 6N.	R	Tetrachloroethylene (Perchloroethylene)	LR			
<b>Bases</b>	Ammonium Hydroxide, 3N.	R	<b>Ketones</b>		Trichloroethylene	LR
	Ammonium Hydroxide, 6N.	R			Acetone	R
	Potassium Hydroxide, 3N.	R		Cyclohexanone	R	
	Sodium Hydroxide, 3N.	R		Methyl Ethyl Ketone	R	
	Sodium Hydroxide, 6N.	R		Methyl Isobutyl Ketone	R	
<b>Alcohols</b>	Amyl	R	<b>Oils</b>	Lubrication Oil	R	
	Benzyl	R		Cottonseed Oil	R	
	Butyl	R		Peanut Oil	R	
	Ethyl	R		Sesame oil	R	
	Isopropyl	R		White Petroleum	R	
	Methyl	R		Lanolin	R	
<b>Esters</b>	Amyl Acetate	R	<b>Organic Solvents</b>	Dimethyl Formamide	R	
	Butyl Acetate	R		Dimethyl Sulfoxide	R	
	Cellulose Acetate	R		Formaldehyde	R	
	Ethyl Acetate	R		Gasoline	LR	
	Isopropyl Acetate	R		Phenol Liquid	R	
	Methyl Acetate	R		Pyridine	R	
<b>Ethers</b>	Diethyl Ether	R	<b>Gases</b>	Turpentine	R	
	Dipropyl Ether	R		Nickel Sulfate	R	
	Dioxane	R		Pentane	R	
	Tetrahydrofuran	R		Helium	R	
				Hydrogen	R	
				Ozone	R	
				Methane	R	
				Nitrogen	R	

**R –Recommended    LR –Limited Recommended    NR –Not Recommended**

This chemical compatibility table is intended for use as a guide only.

Recommendations are based upon static condition of 48 hours at 25°C and 1.0 atmosphere pressure.

## ORDERING GUIDE: KAREI – PTFE – ( A ) – ( B ) – ( C ) – ( D ) – ( E ) – ( F )

<b>( A ) MICRON</b>	005=0.05, 01=0.1, 02=0.2, 04=0.45, 07=0.7, 1=1 um
<b>( B ) TYPE</b>	HB=Hydrophobic membrane, HL=Hydrophilic membrane HLA= ALL PTFE Filter Cartridge
<b>( C ) LENGTH</b>	125, 250, 500, 750, 1000 or 127, 254, 508, 762, 1016 mm
<b>( D ) END STYLE</b>	None=Double Opened End (DOE) S2C=222 & Closed End, S2F=222 & Finned End, S6C=226 & Closed End, S6F=226 & Finned End Note: For SOE with extended adaptor, please include the code of 'EX'.
<b>( E ) GASKET/ O-RING MATERIAL</b>	None=EPDM, V=Viton, S=Silicone, T=Teflon, FEP=Teflon Encapsulated Viton Note : For SOE with stainless steel reinforcement ring, please include the code of 'R'.
<b>( F ) PARTS MATERIAL</b>	None=Polypropylene, RPG=Reinforced PP With Glass, HPE=High Density PE

### EXAMPLE:

- 1) KAREI-PTFE-01HB-250-DOE (PTFE, 0.1 um, Hydrophobic PTFE, 250mm, DOE, EPDM Gasket, P.P. parts material)
- 2) KAREI-PTFE-02HLA-508-S2C (PTFE, 0.2 um, ALL PTFE filter cartridge, 508mm, SOE, 222 Teflon O-Ring, Closed end)
- 3) KAREI-PTFE-01HL-250-S2C-EX-VR-RPG (PTFE, 0.1 um, Hydrophilic PTFE for liquid application, 250mm, SOE, 222 Viton O-Ring with extended adaptor and stainless steel reinforcement ring, Closed end, Reinforced P.P. with glass parts material)

Note: We cannot anticipate all conditions under which this information and our products, or the products of other manufacturers in combination with our products, may be used. We accept no responsibility for results obtained by the applications. Users are advised to make their own testing under actual condition to determine the safety and suitability of each product or product combination for their own purposes and applications. Buyers and users assume all responsibility for liability performance or damage. We reserve the entire right to modify the information without prior notice due to continuous R & D.